




iD Expert Fitting Guide iD Expert Slip

iD Expert Slip

iD Expert Slip offers an effective solution for managing moderate to heavy incontinence, with a high level of absorption and range of sizes.

Ensure the correct size by measuring the user's waist in line with the umbilicus (belly button). Never assume the measurement. Remember, bigger is not better! Then choose the correct absorbency.



- 1 Resealable fixation tapes
- 2 Anti-leak cuffs
- 3 Wetness indicators

Standing

Gently open out the pad and then fold in half lengthways with the absorbent area of the pad facing inwards.

Keeping the pad folded, pass the pad between the legs from front to back.

Open the pad out over the buttocks and lower abdomen, ensuring the side with tapes is placed at the back. Align with the waist and centre correctly using the wetness indicators as a guide.

Ensure the soft elastics are fitted closely into the groin. Fasten the resealable tapes. The position of the tapes is entirely dependent on the user's needs. It is important to create a seal around the legs and groin to reduce risk of leakage.









Do not shake or twist the pad as this can damage it.

Always ensure the tapes are secured to the pad and are not lying on, or close to the skin, as this may cause irritation and discomfort.

When carrying out routine position changes remember to adjust the tapes accordingly to ensure leg seals are maintained.





Lying Down

Gently open out the pad and then fold in half lengthways with the absorbent area of the pad facing inwards.

With the user lying on their side, keep the pad folded and pass between their legs from front to back. Open the pad out over the buttocks and lower abdomen, ensuring the side with the fixation tapes is placed at the back. Align with the waist and centre correctly using the wetness indicators as a guide.

Fasten the resealable tapes on the uppermost side, and then gently roll them and secure the tapes on the opposite side. Now position the user on their back and ensure the anti-leak cuffs are fitted close at the groin. It is important to create a seal around the legs and groin to reduce risk of leakage.









Wetness Indicator

Each pad has two wetness indicator lines to show the presence of urine. These are a guide to when you need to change the product and can be checked without the need to remove the pad.

The yellow line will turn blue and the blue line/writing will fade. When approximately two thirds of the wetness indicator lines have changed colour the pad is ready to be changed.

If there has been any faecal soiling then the pad should be changed immediately.

Important Advice

Oil based barrier creams may reduce the performance of the pads. If barrier creams must be used, they should be applied sparingly and only on treatable areas. Used products should be disposed of in the appropriate manner, they cannot be flushed down the toilet.

Coloured Resealable Tapes (Plastic Backed Only)

There are two main parts to the fixing tapes used on All-In-One products. The blue section is the primary tape and once fixed will remain in place. Attempting to remove this tape once fixed will tear the pad.

The secondary section is the white upper tape which will always remain attached to the back of the pad. This tape can be removed from the blue section and re-fixed to the back of the pad or re-positioned elsewhere on the pad.

This is a simple procedure and involves holding the blue tape with one hand, whilst lifting and removing the white tape with the other hand. You may re-fix the white tape back onto the blue tape as many times as you require.

The Textile back sheet on the iD Expert Slip has single resealable tapes on each side.





Ontex Healthcare UK Weldon House Corbygate Business Park Priors Haw Road Corby NN17 5JG

Phone: 01536 269744

ontexglobal.com

WE LOVE TO HELP YOU CARE

ID EXPERT